
## True 128-bit Extended Validation SSL Certificate

With PaCentric<sup>®</sup>, you can rely on protected and secured communication channels at all times. True 128-bit extended validation. In addition, all patient images are anonymized. PaCentric<sup>®</sup> uses the Internet's most trusted security provider, VeriSign.

The SSL Certificate with extended validation triggers high-security web browsers to display Fimreite Software AS in a green address bar. See image below...



Right click on the PaCentric<sup>®</sup> web site, and then select Properties to get more details.

Current connection: TLS 1.0, AES with 128 bit encryption (High); RSA with 2048 bit exchange.

Properties Certificate		
General		General Details Certification Path
	PaCentric® Web - Medical Image Viewer	Certificate Information
		This certificate is intended for the following purpose(s):
Protocol:	HyperText Transfer Protocol with Privacy	Ensures the identity of a remote computer
Туре:	HTML Document	
Connection:	TLS 1.0, AES with 128 bit encryption (High); RSA with 2048 bit exchange	* Refer to the certification authority's statement for details.
Address: (URL)	https://pacentric.com/pacentric/PaCentric.aspx	Issued to: pacentric.com
Size:	15669 bytes	Issued by: VeriSign Class 3 Extended Validation SSL CA
Created:	17.09.2010	Valid from 05. 02. 2010 to 12. 03. 2011
Modified:	17.09.2010	
	Certificates	Install Certificate Issuer Statement
	OK Cancel Apply	ОК